


Skin and Wound Product Information Sheet

Integra	
Classification	Graft: Dermal Replacement Template
Key Points	<ul style="list-style-type: none"> • Integra dermal regeneration template is a bi-layer membrane (either fenestrated (mesh) or non-fenestrated) used for skin replacement. The base layer (dermal regeneration layer) is composed of bovine tendon collagen and glycosaminoglycan (chondroitin-6-sulfate). The outer layer acts as an epidermal substitute and is made from a thin layer of silicone which helps to control moisture loss from the wound bed and acts as a bacterial barrier. • Application of Integra is a 2-step process: <ul style="list-style-type: none"> ○ It is applied in the operating room (OR) the surgeon places the material over a full thickness surgically excised wound. ○ Over the next 14 - 21 days the cells migrate into the matrix; as the dermal layer matures, the colour of the Integra changes to red, pink, orange, peach and/or cream. The rate of colour change is dependent on the patient’s rate of healing. ○ In approximately 14 - 21 days, the new dermal layer is complete. The silicone layer is then removed in the OR and a skin graft is applied to complete wound closure if needed.
Indications	<ul style="list-style-type: none"> • To provide coverage over a full thickness defect, exposed tendon or bone associated with burn injuries, trauma, cancer or chronic wounds • For the post-excision treatment of full thickness or deep partial thickness injury where sufficient autograft is not available at the time of excision • For reconstruction with scar contracture release
Precautions	<ul style="list-style-type: none"> • The use of Integra in pregnant women has not been clinically evaluated and use should only occur when benefits clearly outweigh risks • Clinically infected wounds must be treated prior to Integra application • Minimize or prevent pressure and shearing forces as these forces can lead to dislodgement of the Integra • Avoid hydrotherapy immersion and protect the dressing when showering. Do not allow water to get under the silicone layer. • Common complications include hematomas, fluid accumulation beneath the silicone layer, infection, premature silicone separation and incomplete take of the Integra.
Contraindications	<ul style="list-style-type: none"> • The use of Integra is contraindicated in patients with known hypersensitivity to bovine collagen or chondroitin material • The bovine components contained in Integra may be a contraindication for application in some cultures • Integra is incompatible with Dakin’s solution, petroleum-based products; e.g. Xeroform, and enzymatic debridement agents
Formats & Sizes	<ul style="list-style-type: none"> • Sheets <ul style="list-style-type: none"> ▪ 5 x 5cm (2x2in) ▪ 10 x 12.5cm (4x5in) ▪ 10 x 25cm (4x10in) ▪ 20 x 25cm (8x10in) <div style="text-align: right;">  </div>
Application Directions	
<p>Integra is applied by the Surgeon in the OR following excision of the wound bed. The Integra can be applied meshed, un-meshed or ‘pie-cruste’d’. It is stapled/sutured and remains in place until the time of autografting (usually 14-21days). There are 4 options for the Integra cover dressing:</p> <ul style="list-style-type: none"> • Non-antimicrobial, non-petroleum interface and NPWT • Bulky dressing and compression • Antimicrobial, bulky dressing and compression • Antimicrobial interface and NPWT 	<p>Rationale</p> <p>Complete excision down to viable tissue is required to ensure Integra takes to the wound bed.</p> <p>The surgeon may use mesh or cut slits in the silicone layer to allow exudate to pass freely from the Integra and prevent fluid accumulation and/or hematoma.</p> <p>Pressure (NPWT or compression) assists with the adherence of the Integra to the underlying wound bed and also minimizes any mechanical dislodgement.</p>



Skin and Wound Product Information Sheet

<p>Care of Integra Layer</p>	
<p>Monitor, and urgently consult with the Surgeon, regarding the following conditions as these will prevent the adherence of Integra to the underlying wound bed. NSWOCs/Wound Clinicians to follow Surgeon’s order for the resolution of the situation:</p> <p><u>Fluid accumulation:</u></p> <ul style="list-style-type: none"> • Evacuate using an 18- or 20-gauge needle and then roll/wick away any remaining fluid with a cotton tipped applicator or gauze pad. <p><u>Hematoma:</u></p> <ul style="list-style-type: none"> • For hematomas in fluid state, evacuate using an 18- or 20-gauge needle and then roll/wick away any remaining fluid with a cotton tipped applicator or gauze pad. • For hematomas no longer in fluid state, Conservative Sharp Wound Debridement is needed, refer to Physician/NP/NSWOC/Wound Clinician. <p><u>Infection:</u></p> <ul style="list-style-type: none"> • Often due to non-viable tissue in the wound bed and/or contamination. The Integra may turn yellow, green, black or brown in color and it may have an odor. • Evacuate purulent drainage using a 18- or 20-gauge needle, incising or rolling the Integra • Treat with a topical antimicrobial dressing: <ul style="list-style-type: none"> ○ Fenestrated(mesh): apply the dressing over the entire sheet. ○ Non-fenestrated: the dressing must be applied to the wound edges, over the Integra seams and any other open areas on the Integra to allow the antimicrobial in under the silicone layer <p><u>Mechanical Shearing:</u></p> <ul style="list-style-type: none"> • Roll or use plastic-coated surfaces for re-positioning the client. If the Integra dressing is on the client’s back, position side-lying or prone (if tolerated), to reduce pressure on the Integra. An air fluidized low loss air bed may be appropriate. • Follow Surgeon’s orders for range-of-motion exercises and ensure the wound is covered securely. 	<p>If using an Acticoat product, Acticoat Flex is to be applied dry; whereas, Acticoat needs to be moistened prior to application.</p>
<p>Outer Dressing Change</p>	
<p>Physician orders are required for NPWT:</p> <ul style="list-style-type: none"> ○ Pressure and therapy settings ○ Type of interface to be used with NPWT ○ Type of wound covering, foam or gauze ○ Dressing change frequency (usually 2-3 times/week) <p><u>Do not</u> pre-soak (instill fluid down the tubing) or use Normal Saline to dampen the dressing for removal.</p> <p>Physician’s orders are required for conventional dressing:</p> <ul style="list-style-type: none"> ○ Type of contact layer, if using ○ Number of gauze dressing/pad layers ○ Type of compression ○ Dressing change frequency (usually 2-3 times/week) <p><u>Do not</u> use Normal Saline to dampen the dressing for removal.</p> <p>Cleansing may cause shearing of the Integra. Cleansing of the wound area is limited to the peri-wound skin and to staples or sutures where crusting has developed.</p>	
<p>Expected Outcome</p>	
<p>New dermal layer will regenerate within 14-21 days. Client will return to the OR for autografting, the 2nd step in the process.</p>	
<p>For further information, please contact your Wound Clinician.</p>	