


Product Information Sheet

Mepore Pro	
Classification	Cover Dressing: Island Dressing
Key Points	<ul style="list-style-type: none"> Comprised of an absorbent wound pad with a non-adherent contact layer, an adhesive border with water based acrylic adhesive, and a breathable and waterproof outer film. Waterproof. Can be used as either a primary or secondary dressing.
Indications	<ul style="list-style-type: none"> For low to moderately exuding wounds.
Precautions	<ul style="list-style-type: none"> Do not use on fragile skin as may cause skin irritation. Has not been evaluated on pregnant/lactating individuals or neonates/infants, consult with physician/NP prior to using on these populations.
Contraindications	<ul style="list-style-type: none"> Sensitivity or allergy to acrylic adhesive or other components of the dressing.
Formats & Sizes	<ul style="list-style-type: none"> Adhesive pad <ul style="list-style-type: none"> 6 x 7 cm 9 x 10 cm 9 x 15 cm 9 x 20 cm 9 x 25 cm 9 x 30 cm



Directions	Rationale / Key Points
Selection	
<p>Select appropriately sized Mepore Pro dressing.</p> <p>For a wound with depth, choose appropriate wound filler for amount of exudate expected and anticipated frequency of dressing change.</p>	Ensure absorbent pad covers the wound.
Preparation	
<p>Cleanse wound and periwound/ surrounding skin with sterile normal saline or agency approved wound cleanser.</p> <p>Dry periwound / surrounding skin.</p> <p>If required, apply barrier film to periwound skin. Allow to thoroughly dry before applying dressing.</p>	<p>See Wound Cleansing Procedure or QR Code below.</p> <p>Barrier film may be used to protect periwound skin from moisture associated skin damage (MASD) and medical adhesive related skin injury (MARSI). Do not use moisturizers, gels or ointments on periwound /surrounding skin as these may interfere with adhesion.</p>
Application	
<p>For wounds with minimal depth (less than 1 cm): remove the release liners and apply sticky side to the skin without stretching. Ensure that absorbent pad covers wound.</p> <p>For wounds with depth (more than 1 cm): lightly fill the dead space up to skin level with appropriate wound filler. Remove the release liners and apply sticky side to the skin without stretching. Ensure that absorbent pad covers wound.</p> <p>Gently press border to ensure proper adhesion to skin.</p>	<p>To protect periwound / surrounding skin from medical adhesive related skin injury (MARSI):</p> <ul style="list-style-type: none"> Do not stretch dressing when applying, Do not apply under tension. Apply over joints in such a way as to not impair flexibility. Avoid using when swelling is expected in the area.
Removal	
<p>Consider using adhesive remover to assist with removal.</p> <p>Grasp edge of dressing and slowly peel dressing from skin in direction of hair growth.</p>	<p>To decrease risk of medical adhesive related skin injury.</p> <p>To minimize discomfort with removal.</p>

Product Information Sheet

Directions	Rationale / Key Points
Frequency of Dressing Change	
Change before dressing is fully saturated. May be left in place for up to 7 days.	Dressing change frequency is dependent on amount of wound exudate.
Expected Outcomes	
Exudate is managed with no periwound maceration. Product performs as expected.	If product does not perform as expected, notify NSWOC/Wound Clinician and then consider submitting a Supply Chain Product Concern Form .
QR Codes	
	
Wound Cleansing Procedure	
For further information please contact NSWOC/Wound Clinician	