



Skin and Wound Product Information Sheet

Prisma	
Classification	Biochemical Modifier with Silver: Inhibitor- Bioabsorbable
Key Points	<ul style="list-style-type: none"> • Prisma Wound Balancing Matrix contains 55% bovine collagen, 44% plant-based Oxidized Regenerated Cellulose (ORC), and 1% ionic silver. • Inactivates excess proteases, binds and protects endogenous growth factors, reduces levels of pro-inflammatory cytokines, scavenges free radicals, enhances cell proliferation. • Product is bio-absorbable. • Product has antimicrobial action.
Indications	<ul style="list-style-type: none"> • To correct the cellular imbalance in a ‘stalled’ healable wound with a clean wound bed that show signs and symptoms (S&S) of local wound infection.
Precautions	<ul style="list-style-type: none"> • Consult with NSWOC/Wound Clinician prior to using this product. • May be used in conjunction with compression therapy.
Contraindications	<ul style="list-style-type: none"> • Do not use for clients with known sensitivity or allergy to cellulose, collagen or silver. • Do not use in the presence of necrotic tissue. • Do not use on 3rd degree burns.
Formats & Sizes	<ul style="list-style-type: none"> • Sheet <ul style="list-style-type: none"> ▪ 28 cm² ▪ 123 cm² <div style="text-align: right;">  </div>
Application Directions	
Cleanse/irrigate wound with sterile normal saline or agency approved wound cleanser; dry peri-wound skin.	Reduces wound debris and allows for adhesion of dressing or tape.
To Apply	
Cut Prisma to completely cover the wound bed.	
If the wound is moist, apply dry piece of Prisma directly on the wound.	When Prisma absorbs exudate or is moistened, it forms a gel. Prisma needs moisture to be effective.
If the wound is dry, apply Prisma on the wound and then moisten it lightly with normal saline.	
Apply appropriate cover dressing to maintain a moisture-balanced wound environment.	The choice of cover dressing is depended upon the amount of exudate expected.
To Remove	
Dressing will have turned into a gel. Prisma is biodegradable, there will be no residual product to remove.	If product is present at dressing change then increase the time between dressing changes to allow it to fully absorb.
Frequency of Dressing Change	
Will depend upon the amount of exudate. Prisma will usually be absorbed over 3 days.	Prisma may need replacing more frequently if the wound is heavily exuding.
Expected Outcome	
S&S of local wound infection are resolved within 2 weeks.	
Measurable improvement in wound size within 3 weeks.	
For further information, please contact your NSWOC/Wound Clinician.	