




## Skin and Wound Product Information Sheet

<b>Promogran</b>	
<b>Classification</b>	<b>Biochemical Modifier: Inhibitor – Bioabsorbable</b>
<b>Key Points</b>	<ul style="list-style-type: none"> <li>• Promogran contains 55% bovine collagen and 45% plant-based Oxidized Regenerated Cellulose (ORC).</li> <li>• Inactivates excess proteases, binds and protects endogenous growth factors, reduces levels of pro-inflammatory cytokines, scavenges free radicals and enhances cell proliferation.</li> <li>• Product is bio-absorbable.</li> <li>• Product has hemostatic properties.</li> </ul>
<b>Indications</b>	<ul style="list-style-type: none"> <li>• To correct the cellular imbalance in ‘stalled’, clean, healable wounds.</li> <li>• May be used under compression therapy.</li> </ul>
<b>Precautions</b>	<ul style="list-style-type: none"> <li>• Consult with Wound Clinician prior to using this product.</li> </ul>
<b>Contraindications</b>	<ul style="list-style-type: none"> <li>• Do not use for clients with known sensitivity or allergy to cellulose or collagen.</li> <li>• Do not use in the presence of necrotic tissue, heavy colonization or infection.</li> </ul>
<b>Formats &amp; Sizes</b>	<ul style="list-style-type: none"> <li>• Sheet                             <ul style="list-style-type: none"> <li>▪ 28 cm<sup>2</sup></li> <li>▪ 123 cm<sup>2</sup></li> </ul> </li> </ul> <div style="text-align: right;">  </div>
<b>Application Directions</b>	
Cleanse/irrigate wound with sterile normal saline or agency approved wound cleanser, dry peri-wound skin.	Reduces wound debris and allows for adhesion of dressing or tape.
<b>To Apply</b>	
Cut Promogran to completely cover the wound bed.	
If the wound is moist, apply dry piece of Promogran directly on the wound.	When Promogran is moistened or absorbs exudate, it forms a gel.
If the wound is dry, apply Promogran on the wound and then moisten it lightly with normal saline.	Promogran needs moisture to be effective.
Cover with appropriate moisture retaining cover dressing.	The choice of cover dressing is depended upon the amount of exudate expected.
<b>To Remove</b>	
Promogran is biodegradable, there will be no residual product to remove.	If product is present at dressing change then increase the time between dressing changes to allow it to fully absorb.
<b>Frequency of Dressing Change</b>	
Will depend upon to amount of exudates. Promogran will usually be absorbed over 3 days.	Promogran may need replacing more frequently if the wound is heavily exuding.
<b>Expected Outcome</b>	
Measurable improvement in wound size within 3 weeks.	
<b>For further information, please contact your NSWOC/Wound Clinician.</b>	