



Provincial Professional Practice Stream Wound Ostomy Continence Document Update Notification of New, Revised, Deleted – Master			
Aug 21, 2024	Product Information Sheet	<a href="#">OASIS Bedbath Wipes</a>	New document.
Aug 19, 2024	Product Information Sheet	<a href="#">Melgisorb Ag</a>	Revised: <ul style="list-style-type: none"> <li>• Added consult NSWOC/WC prior to using product in deep undermining/tunneling.</li> <li>• Added: product shrinks as it absorbs exudate so a larger sheet or more than one sheet may be needed</li> <li>• Added the 3cm wide ribbon should not be cut lengthwise, (e.g., make a 1.5cm wide ribbon) due to shrinkage &amp; compromised tensile strength.</li> <li>• Added for wounds with heavy exudate daily or every second day dressing changes may be needed.</li> </ul> Document reformatted.
Aug 18, 2024	Product Information Sheet	<a href="#">Anasept Skin &amp; Wound Cleanser</a>	Revised: <ul style="list-style-type: none"> <li>• May be used in conjunction with maggot therapy.</li> <li>• Removed the need for doing a normal saline/sterile water rinse prior to applying an antimicrobial dressing; wound bed just needs to be pat dry to remove excess solution.</li> </ul> Document reformatted.
Aug 18, 2024	Product Information Sheet	<a href="#">Vashe Wound Solution</a>	Revised/reformatted. As per vendor: <ul style="list-style-type: none"> <li>• Removed the need for doing a normal saline/sterile water rinse prior to applying an antimicrobial dressing, wound bed just needs to be pat dry to remove excess solution.</li> <li>• Revised solution stability to now be until the bottle (opened or unopened) reaches is shelf-life date.</li> </ul> Document reformatted.
Aug 14, 2024	Product Information Sheet	<a href="#">Medigrip LF (Lower Limb Application)</a>	Revised/updated: <ul style="list-style-type: none"> <li>• Revised the BC Practice statement: now calls out that for compression less than 20mmHg ABPI and/or TPBI are not required given pedal pulses are palpable and cap refill is normal.</li> <li>• Updated sizing chart to just show measurements in centimeters.</li> </ul> Document reformatted.



Aug 14, 2024	Poster	<a href="#">Medigrip LF (Lower Limb Application Sizing Chart)</a>	New: chart that lays out what size of Medigrip LF bandage to use based upon calf measurement and compression required.
Aug 5, 2024	Product Information Sheet	<a href="#">Dakin's Solution</a>	Revised: change made to call out using only non-concentrated household/laundry sodium hypochlorite-based bleach, such as Clorox or Javex, to ensure a 5.25% concentration.
Jul 11, 2024	Quick Reference Guide	<a href="#">Wound Management for Frostbite Injury: Guide</a>	Revised: title changed to call out that the guide speaks only to the wound management of a frostbite injury.
Jul 11, 2024	Guideline	<a href="#">Wound Management for Frostbite Injury: Guideline</a>	Revised: title changed to call out that the guideline speaks only to the wound management of a frostbite injury.
Jun 24, 2024	Guideline/Procedure	<a href="#">High Voltage Pulsed Current Wound Therapy: Guideline/ Procedure</a>	Updated: Silverlon dressing now back in Canada and on the WOC Formulary. Document updated re the use of Silverlon dressings as a HVPC conductor.
Jun 21, 2024	Clinical Resource	<a href="#">Yukon Ostomy Suppliers List</a>	Updated.
Jun 21, 2024	Product Information Sheet	Biatain Ag Adhesive	Removed: vendor discontinuation.
Jun 14, 2024	Product Information Sheet	<a href="#">Acticoat Flex 3/Acticoat Flex 7</a>	Updated/reformatted.
Jun 14, 2024	Product Information Sheet	<a href="#">Acticoat/Acticoat 7</a>	Updated/reformatted.
Jun 14, 2024	Product Information Sheet	Acticoat Flex 3/Acticoat Flex 7 for Burns	Removed: Updated PISheet now for both wounds & burns.
Jun 14, 2024	Product Information Sheet	Acticoat Flex 3/Acticoat Flex 7 for Wounds	
Jun 14, 2024	Product Information Sheet	Acticoat/Acticoat 7 for Burns	Removed: Updated PISheet now for both wounds & burns.
Jun 14, 2024	Product Information Sheet	Acticoat/Acticoat 7 for Wounds	
Jun 14, 2024	Product Information Sheet	<a href="#">CURAD Non-adherent Oil Emulsion Dressing</a>	New: backorder sub for the Adaptic roll (7.6x152.4cm)
Jun 14, 2024	Product Information Sheet	Ag+Powder	Removed: product discontinued by vendor.
May 28, 2024	Additional Education & Competencies	<a href="#">NPWT Level Two: Dressing Application Learning Plan</a>	Updated: Added there are six videos to view (approx. 25mins for all) once the module has been completed.
May 14, 2024	Video	<a href="#">Setting Up an ActiVAC Device</a>	New: six short videos that support the NPWT Level Two: Open Wound Dressing Application module.
May 14, 2024	Video	<a href="#">Setting Up a VACUIta Device</a>	



May 14, 2024	Video	<a href="#">Y-Connecting Two Wounds</a>	
May 14, 2024	Video	<a href="#">Offloading the Trac-Pad</a>	
May 14, 2024	Video	<a href="#">Bridging Two or More Wounds</a>	
May 14, 2024	Video	<a href="#">Dressing Application for VACUta &amp; ActiVAC Devices</a>	
May 04, 2024	Clinical Resource	<a href="#">BC Ostomy Suppliers – Vancouver Island</a>	Updated.
May 04, 2024	Clinical Resource	<a href="#">BC Ostomy Suppliers – Northern BC</a>	Updated.
May 04, 2024	Clinical Resource	<a href="#">BC Ostomy Suppliers – Lower Mainland/ Sunshine Coast</a>	Updated.
May 04, 2024	Clinical Resource	<a href="#">BC Ostomy Suppliers – Interior BC</a>	Updated.
May 04, 2024	Clinical Resource	<a href="#">BC Ostomy Suppliers – Benefits and Supplier Info</a>	Update with May date.
May 02, 2024	Clinician Resource	<a href="#">Mattress Audit Process</a>	New document.
Mat 02, 2024	Clinician Resource	<a href="#">Mattress Audit Tool</a>	New document.
Apr 19, 2024	Quick Reference Guide	<a href="#">Recommended LTC S&amp;W Product List</a>	Updated.
Mar 07, 2024	Product Information Sheet	<a href="#">Jelonet</a>	Updated/reformatted.
Mar 06, 2024	Product Information Sheet	<a href="#">Bactigras</a>	Updated/reformatted.
Feb 29, 2024	Product Information Sheet	<a href="#">Cavilon No Sting Barrier Film</a>	Revised: added in wording re using product for doing a crusting and added a 'how to do' crusting technique section.
Feb 29, 2024	Product Information Sheet	<a href="#">Essential Wipes</a>	Revised:
Feb 29, 2024	Product Information Sheet	<a href="#">Comfort Wipes</a>	Revised: identified Demethicone % as 3%.
Feb 21, 2024	Clinician Resource	<a href="#">BC Contenance Clinics</a>	Updated: direct phone line for Abbotsford clinic.
Feb 20, 2024	Product Information Sheet	<a href="#">AtmosAir 9000 Series</a>	Updated/reformatted: AtmosAir 9000AR (lateral rotation) surface is no longer manufactured.
Feb 19, 2024	Product Information Sheet	<a href="#">Ag+ Powder</a>	Updated/reformatted: merged Ag+ Powder for Skin and Ag+ Powder for Wounds PISheets into one PISheet
Feb 19, 2024	Product Information Sheet	<a href="#">Inadine</a>	Updated/reformatted.
Feb 19, 2024	Product Information Sheet	<a href="#">Biatin Ibu</a>	Updated/reformatted.



Feb 19, 2024	Product Information Sheet	<a href="#">No Sting Skin Prep</a>	Re-posted.
Feb 19, 2024	Product Information Sheet	Turning Clock	Removed: no longer being manufactured.
Feb 02, 2024	Client Health Education Resource	<a href="#">Taking Urostomy Urine Specimen for Culture &amp; Susceptibility (C&amp;S)</a>	New document.
Feb 02, 2024	Clinical Resource	<a href="#">BC Health Authority Outpatient/Ambulatory Clinics – Total Contact Casting (TCC)</a>	Updated: added VCH Raven Song Community Clinic to the list.
Jan 19 <sup>th</sup> , 2024	Quick Reference Guide	<a href="#">Skin Care Protocol – Adults</a>	Updated: add a photo of DriGo-HP.
Jan 15 <sup>th</sup> , 2024	Product Information Sheet	<a href="#">Deodorant Clean-Up Washcloths</a>	Updated: removed wording re dimethicone as product doesn't contain dimethicone.
Jan 15 <sup>th</sup> , 2024	Product Information Sheet	CalciCare Ag	Removed: no longer on contract.
Nov 29 <sup>th</sup> , 2023	Decision Support Tool	<a href="#">Frostbite Injury: Quick Reference Guide</a>	Updated: <ul style="list-style-type: none"> <li>Added new aloe-vera product Aloex to the documents as the Aloe Vesta product has been discontinued by the vendor</li> <li>The aloe-vera product can be used beyond the first five days of the post-injury period.</li> </ul>
Nov 29 <sup>th</sup> , 2023	Decision Support Tool	<a href="#">Frostbite Injury: Guideline</a>	
Nov 27 <sup>th</sup> , 2023	Product Information Sheet	<a href="#">Aloex Pure Gel for Frostbite Treatment</a>	New document.
Nov 27 <sup>th</sup> , 2023	Product Information Sheet	Aloe Vesta Protective Ointment for Frostbite Treatment	Removed, vendor discontinuation.
Oct 30, 2023	Assessment Form	<a href="#">Community Care Established Ostomy Assessment Flow Sheet &amp; Management Plan</a>	New document.
Oct 30, 2023	Documentation Guide	<a href="#">Community Care Established Ostomy: Documentation Guide</a>	New document.
Oct 26, 2023	Product Information Sheet	<a href="#">Melgisorb Ag</a>	New document.
Oct 17, 2023	Quick Reference Guide	<a href="#">Skin Care Protocol - Adults</a>	Revised: <ul style="list-style-type: none"> <li>Renamed Remedy skin care products</li> <li>Replaced Interdry-Ag with DriGo-HP.</li> </ul>
Oct 11, 2023	Product Information Sheet	<a href="#">DriGo-HP</a>	New document.
Oct 11, 2023	Product Information Sheet	<a href="#">Cavilon No Sting Barrier Film</a>	New document.



Oct 11, 2023	Product Information Sheet	No Sting Skin Prep	Removed: product no longer on contract.
Oct 06, 2023	Decision Support Tool	<a href="#">Nutrition for Wound Prevention &amp; Management: Guideline for Dietitians</a>	New document.
Oct 06, 2023	Decision Support Tool	<a href="#">Negative Pressure Wound Therapy: Guideline</a>	Revised: <ul style="list-style-type: none"> <li>• ‘Reusable’ changed to ‘multi-use’.</li> <li>• ‘Disposable’ changed to ‘single-use’.</li> <li>• Clarified that HA policy needs to be in place for NSWOCs/WC to provide an order for NPWT.</li> <li>• Clarified guidance re apply/reapplying Closed Incision NPWT dressings.</li> <li>• Added NPWTi-d Smart Instill info to the NPWT order section.</li> <li>• Clarified checking NPWT orders for the Therapy Settings only needs to be done at the beginning of each shift and that documentation of the q2h safety/ monitoring check is as per HA policy.</li> </ul>
Oct 06, 2023	Module & Videos	<a href="#">NPWT Level Two: Dressing Application – Closed Incision &amp; Skin Graft Dressing Management</a>	New module
Oct 06, 2023	Module & Videos	<a href="#">NPWT Level Two: Dressing Application – Open Wound Principles</a>	New module.
Oct 06, 2023	Education Requirements & Competencies	<a href="#">NPWT Level Two: Dressing Application Learning Plan</a>	New document.
Oct 06, 2023	Education Requirements & Competencies	<a href="#">NPWT Level Two: Dressing Application Competencies</a>	New document.
Sep 24, 2023	Decision Support Tool	<a href="#">NPWT Dressing Application for VACUta4 &amp; ActiVAC: Procedure</a>	Revised dressing procedure: <ul style="list-style-type: none"> <li>• Changed ‘reusable’ to ‘multi-use’</li> <li>• Renamed ‘window-paning’ to ‘picture-framing’ as per learning module.</li> <li>• Aligned the measurements used in the procedure and the learning module re length and wide of peri-incisional protection dressing, interface, black foam and transparent drape.</li> </ul> Revised Y-Connector procedure: <ul style="list-style-type: none"> <li>• Added in statement to use the male end of the Y-connector for the wound of most concern.</li> <li>• Not to connect wounds of different etiologies</li> </ul>



			<ul style="list-style-type: none"> <li>• Not to be used with Veraflo.</li> </ul>
Sep 24, 2023	Decision Support Tool	<a href="#">Veraflo NPWT Dressing Application: Procedure</a>	Revised: <ul style="list-style-type: none"> <li>• New Smart Instill steps added to document.</li> <li>• As per vendor, cassette now to be changed q3days and whenever the type of solution was changed; previously was every 11,250ml of solution used or every 7 days, whichever comes first and whenever the type of solution is changed.</li> </ul>
Sep 21, 2023	Quick Reference Guide	<a href="#">Adult Skin Care Protocol</a>	Updated: revised with new names and photos for Remedy Cleanse, Moisturize & Prevent products.
Sep 18, 2023	Product Information Sheet	Secura Dimethicone Protectant	Removed: old resource.
Sep 18, 2023	Product Information Sheet	Secura Professional Care Lotion	Removed: old resource.
Sep 18, 2023	Product Information Sheet	Secura Foam Cleanser	Removed: old resource.
Sep 18, 2023	Product Information Sheet	<a href="#">Remedy Cleanse Hydrating Shampoo &amp; Body Wash</a>	New resource.
Sep 18, 2023	Product Information Sheet	Remedy Shampoo & Body Wash	Removed: name change.
Sep 18, 2023	Product Information Sheet	<a href="#">Remedy Prevent</a>	New resource.
Sep 18, 2023	Product Information Sheet	Remedy Hydraguard	Removed: name change.
Sep 18, 2023	Product Information Sheet	<a href="#">Remedy Cleanser</a>	New resource.
Sep 18, 2023	Product Information Sheet	Remedy Foam Cleanser	Removed: name change.
Sep 13, 2023	Product Information Sheet	<a href="#">Remedy Moisturize</a>	New resource.
Sep 13, 2023	Product Information Sheet	Remedy Nourishing Cream	Removed: name change.
Jul 06, 2023	Clinical Resource	<a href="#">QuickChange Wrap Application: Edematous</a>	New resource.
Jul 06, 2023	Clinical Resource	<a href="#">QuickChange Wrap Application: Hidden</a>	New resource.
Jul 06, 2023	Clinical Resource	<a href="#">QuickChange Wrap Application: Short/Retracted</a>	New resource.
Jul 06, 2023	Clinical Resource	<a href="#">QuickChange Wrap Application: Standard/Long</a>	New resource.
Jul 06, 2023	Product Information Sheet	<a href="#">QuickChange Wrap</a>	New resource.



Jul 06, 2023	Product Information Sheet	<a href="#">Mepilex Border Flex Lite</a>	New resource.
Jul 06, 2023	Product Information Sheet	Mepilex Border Lite	Removed: product discontinued.
Jul 03, 2023	Modules & Videos	Semiology of Cutaneous Lesions: Dermatology Education Session#1	Removed.
Jul 03, 2023	Modules & Videos	Diagnosis of a Cutaneous Injury/Ulceration & Creation of an Action Plan - Dermatology Education Session#2	Removed.
June 9 <sup>th</sup> , 2023	Videos	<a href="#">NPWT Video#1 Safety Monitoring Check</a> <a href="#">NPWT Video#2 Trouble Shooting an Air Leak</a> <a href="#">NPWT Video#3 Therapy Discontinued/Completed: Removal of Dressing – Closed Incision&amp; Skin Graft</a> <a href="#">NPWT Video#4 Therapy Discontinued/Completed: Removal of Dressing – Open Wound</a> <a href="#">NPWT Video#6: Irreparable Air Leak Device Failure: Alternate Dressing -Take Down</a> <a href="#">NPWT Video#5 Irreparable Air Leak Device Failure: Alternate Dressing – Cut &amp; Cover</a>	Revised: each video reviewed and video tightened up to shorten its length where possible; new voice-over done.
June 9th, 2023	Module	NPWT Level One: Monitor & Manage Module	Revised: ensured not to repeat information that is in the Learning Plan. Tighten-up wording for the Adverse Events screen as well as added in what to do in the event of S&S of wound infection. Reviewed/ revised Learning Activities to strengthen the learning of what to do in the event of an alert/alarm. These learning activities and the knowledge quiz should prepare the nurse to do a safety/monitoring check without the need for a return demo with Mentor (unless needed as per Mentor).



May 20, 2023	Education Requirements & Competencies	<a href="#">NPWT Level One: Monitor &amp; Manage Competencies</a>	Revised: reformatted to set up return demonstration of safety/monitoring check to be at the CNE/Nurse Mentor discretion. Moved previous Skill statements re identifying where to find information on alerts/alarms and canister change to Knowledge section as well as adding statement re identifying where to find information on NPWT device safety considerations.
May 20, 2023	Education Requirements & Competencies	<a href="#">NPWT Level One: Monitor &amp; Manage Lesson Plan</a>	Revised: removed the requirement of doing a return demonstration for how to do and document a safety/monitoring check.
May 19, 2023	Mini - Product Information Sheet	<a href="#">PolyMem</a>	New document: backorder sub for two sizes of Adaptic Digit.
May 11, 2023	Product Information Sheet	<a href="#">EdemaWear</a>	Revised: in the Key Points/Compression Therapy section, removed all references to compression wraps as EdemaWear is a stocking.
May 11, 2023	Product Information Sheet	<a href="#">Tubigrip (Lower Leg Application)</a>	Revised: in the Key Points/Compression Therapy section, removed all references to compression wraps as Tubigrip is a stocking.
May 11, 2023	Decision Support Tool Procedure	<a href="#">Urostomy Urine Collection for C&amp;S: Procedure</a>	New document.
May 11, 2023	Clinical Resource	<a href="#">Ostomy Patient Peer Support Information Sheet</a>	New document.
May 01, 2023	Clinical Resource	<a href="#">BC Sepsis Network Resources</a>	New: link to the BC Sepsis Network which has clinician resources for identifying/ treating sepsis and post-sepsis syndrome as well as patient/family resources for dealing with post-sepsis syndrome. The link 'lives' in the Wound Infection folder.
Apr 11, 2023	Modules & Videos	<a href="#">Diagnosis of a Cutaneous Injury/Ulceration &amp; Creation of an Action Plan - Dermatology Education Session#2</a>	New video: recording of the April 6 <sup>th</sup> Dermatology Education Session#2.
Apr 03, 2023	Decision Support Tool	<a href="#">LPN: Care of a Wound with a Non-Visible Wound Bed AER&amp;C</a>	Revised: corrected the link to the WATFS Documentation Guide as well as removed the quiz from the appendix as quiz is on LHub.
Apr 03, 2023	Decision Support Tool	<a href="#">Skin Care Protocol</a>	Revised: added in the use of Remedy Foam Cleanser as a perineal cleanser.





Mar 27, 2023	Modules & Videos	<a href="#">Semiology of Cutaneous Lesions: Dermatology Education Session#1</a>	New video: recording of the March 9 <sup>th</sup> Dermatology Education Session#1.
Mar 23, 2023	Assessment Form	<a href="#">Braden Risk &amp; Skin Assessment Flow Sheet</a>	Revised: with permission from Dr. Nancy Bergstrom, the Braden Scale for Predicting Pressure Sore Risk form's Nutrition sub-scale was revised to have 'protein-rich foods' wording with meat and dairy as examples (previous wording was just meat and dairy).
Mar 23, 2023	Decision Support Tool	<a href="#">Frostbite Injury: Guideline</a>	Revised: updated the Practice Level to clarify the LPN scope.
Mar 23, 2023	Decision Support Tool	<a href="#">Hypergranulation Tissue: Guideline</a>	Revised: updated the Practice Level to clarify the LPN scope.
Mar 20, 2023	Mini - Product Information Sheet	<a href="#">Unna-Z</a>	New document: backorder sub for Viscopaste.
Feb 28, 2023	Clinical Resource	<a href="#">BC PharmaCare Enrolled Orthotics Providers: List</a>	New document.
Feb 28, 2023	Clinical Resource	<a href="#">BC Health Authority Outpatient/Ambulatory TCC Clinics: List</a>	New document.
Feb 28, 2023	Clinical Resource	<a href="#">Orthotics Referral Post Completion TCC: Form</a>	New document.
Feb 28, 2023	Clinical Resource	<a href="#">BC PharmaCare Coverage TCC Offloading Devices: Process</a>	New document.
Feb 21, 2023	Assessment Form	<a href="#">Frostbite Injury Assessment Form</a>	Revised: document now supports up to six assessments.
Feb 21, 2023	Decision Support Tool	<a href="#">Frostbite Injury: Quick Reference Guide</a>	Revised: as per guideline.
Feb 21, 2023	Decision Support Tool	<a href="#">Frostbite Injury: Guideline</a>	Revised: last of the Provincial Partners Review feedback incorporated.
Feb 21, 2023	Product Information Sheet	<a href="#">Tubigrip (Lower Leg Application) Sizing Chart</a>	New: poster of the Tubigrip sizing chart.
Feb 21, 2023	Product Information Sheet	<a href="#">Tubigrip (Lower Leg Application)</a>	Revised: sizing chart revised to only display measurements in centimeter. Moved statement re PISheet not discussing application for sprains, post-burn compression to the top of the Key Point section (was in the Indications section previously).
Jan 28, 2023	Product Information Sheet	<a href="#">ABENA Pants</a>	Revision: product renamed from Abri-Flex; hip measurements are to be used to determine size of the pants (not waist and hip). PISheet now shows measurements only in centimeters.

British Columbia Provincial Professional Practice Stream – Wound Ostomy Continence



Jan 28, 2023	Product Information Sheet	<a href="#">ABENA San</a>	Revision: product renamed from Abri-San Protective Pad; small adjustments made to the length and width of the pads. PISheet now shows measurements only in centimeters.
Jan 18, 2023	Product Information Sheet	<a href="#">TotalCare Bariatric Plus Therapy System</a>	New document.
Jan 18, 2023	Decision Support Tool - Guideline	<a href="#">Hypergranulation Tissue: Guideline</a>	Revised: in the DST Indications for Use section, the term 'Vascular Access Device' changed to 'Ventricular Assisted Device'.
Jan 13, 2023	Product Information Sheet	<a href="#">Jelonet</a>	Updated: 10cm x 7m dressing & 10 x 10cm tray discontinued.
Jan 11, 2023	Product Information Sheet	<a href="#">Therapulse ATP</a>	New document.